

AGROBACTERIUM COMPETENT CELLS

Day 1

AGROBACTERIUM COMPETENT CELLS

Solutions: * Autoclave solutions 35min

NaCl 0.15M 500ml = 15mL NaCl 5M + 485 ml H₂O

CaCl₂ 20mM 200ml = 5.88mg CaCl₂ + 200ml H₂O

Keep solutions at 4°C

LB liquid media

Protocol

Day 1: Streak agrobacterium strain GV3101 from a -80°C aliquot into a LB +Rif (100ug/ml) plate

Day 2 (night): Prepare a 50ml liquid culture (LB+Rif) from the plate and grow it overnight

Day 3 (early morning): 5 ml from the saturated o/n culture into a 100ml LB+Rif media and let it grow to a OD =1.65

Day 3 (evening): ICE ALL THE TIME

- Pellet the cells by using a centrifuge: 5min 4°C 6000rpm
 - Wash pellet with 30ml of sterile NaCl 0.15M
 - Pellet the cells by using a centrifuge: 5min 4°C 6000rpm
 - Resuspend pellet with 2mL CaCl₂ 20mM
 - make 100ul aliquots and freeze with N₂